

XLOC[®] POLY RADIUS

PLATES & SCREWS



TABLE OF CONTENTS

Table of Contents	3
XLOC® Screws	4
Radius XLOC® Poly plates	6
Radius XLOC® Poly Plates Set	9
Radius XLOC® Poly Instrument Set	10

XLOC® SCREWS

XLOC® screws are made of high-quality, Type II anodized Titanium and offer excellent biocompatibility, durability and corrosion resistance. They are de-

signed for orthopedic procedures and ensure stable fixation. The T10 screw head ensures safe handling and controlled force transmission.

PRODUCT FEATURES


Type II Anodized

Enhanced durability and reliability

The increased surface hardness significantly improves fatigue strength, effectively preventing wear and reducing the risk of cold welding.

Superior protection and compatibility

Excellent corrosion resistance combined with proven biocompatibility ensures long-lasting performance and safe use.



Cortical screws

2.7 mm | T10 | Self-cutting | Titanium



ITEM NO.	LENGTH
16-702710	10 mm
16-702712	12 mm
16-702714	14 mm
16-702716	16 mm
16-702718	18 mm
16-702720	20 mm
16-702722	22 mm
16-702724	24 mm
16-702726	26 mm
16-702728	28 mm
16-702730	30 mm
16-702732	32 mm
16-702734	34 mm
16-702736	36 mm
16-702738	38 mm
16-702740	40 mm
16-702742	42 mm
16-702744	44 mm
16-702746	46 mm
16-702748	48 mm
16-702750	50 mm

Special Feature

Reliable and versatile; the tip enables self-tapping, precise, and time-efficient insertion

Application

Stable and long-lasting fixation

XLOC® Screws

2.7 mm | T10 | Self-cutting | Titanium



ITEM NO.	LENGTH
16-812710	10 mm
16-812712	12 mm
16-812714	14 mm
16-812716	16 mm
16-812718	18 mm
16-812720	20 mm
16-812722	22 mm
16-812724	24 mm
16-812726	26 mm
16-812728	28 mm
16-812730	30 mm
16-812732	32 mm
16-812734	34 mm
16-812736	36 mm
16-812738	38 mm
16-812740	40 mm
16-812742	42 mm
16-812744	44 mm
16-812746	46 mm
16-812748	48 mm
16-812750	50 mm

Special Feature

Reliable and angle-stable. The locking mechanism ensures consistent stability.

Application

Precise and reliable fixation in orthopedic procedures

XLOC® Poly Screws

2.7 mm | T10 | Self-cutting | Titanium



ITEM NO.	LENGTH
16-822710	10 mm
16-822712	12 mm
16-822714	14 mm
16-822716	16 mm
16-822718	18 mm
16-822720	20 mm
16-822722	22 mm
16-822724	24 mm
16-822726	26 mm
16-822728	28 mm
16-822730	30 mm
16-822732	32 mm
16-822734	34 mm
16-822736	36 mm
16-822738	38 mm
16-822740	40 mm
16-822742	42 mm
16-822744	44 mm
16-822746	46 mm
16-822748	48 mm
16-822750	50 mm

Special Feature

Angulation up to 15° in all directions

Application

Adaptable to individual fracture patterns

RADIUS XLOC® POLY PLATES

The Radius XLOC® Poly bone plates (Type II anodized) combine maximum stability with minimal weight. The biocompatible material ensures excellent tolerability and promotes optimal bone healing. Thanks to their anatomically pre-shaped design, the plates

naturally conform to curved bone structures and enable precise and secure fixation. Ideal for demanding surgical applications where reliability and long-term performance are critical.

PRODUCT FEATURES


Type II Anodized

Enhanced durability and reliability

The increased surface hardness significantly improves fatigue strength, effectively preventing wear and reducing the risk of cold welding.

Superior protection and compatibility

Excellent corrosion resistance combined with proven biocompatibility ensures long-lasting performance and safe use.



Left

Radius XLOC® Poly Plates

Distal | Head width 19.5 mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80162	6	2	45 mm	19.5 mm
11-80163	6	3	53 mm	19.5 mm
11-80164	6	4	61 mm	19.5 mm
11-80165	6	5	69 mm	19.5 mm
11-80166	6	6	77 mm	19.5 mm
11-80167	6	7	85 mm	19.5 mm
11-80168	6	8	93 mm	19.5 mm
11-80169	6	9	101 mm	19.5 mm
11-80170	6	10	109 mm	19.5 mm
11-80171	6	11	117 mm	19.5 mm

Right

Radius XLOC® Poly Plates

Distal | Head width 19.5 mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80062	6	2	45 mm	19.5 mm
11-80063	6	3	53 mm	19.5 mm
11-80064	6	4	61 mm	19.5 mm
11-80065	6	5	69 mm	19.5 mm
11-80066	6	6	77 mm	19.5 mm
11-80067	6	7	85 mm	19.5 mm
11-80068	6	8	93 mm	19.5 mm
11-80069	6	9	101 mm	19.5 mm
11-80070	6	10	109 mm	19.5 mm
11-80071	6	11	117 mm	19.5 mm

Left

Radius XLOC® Poly Plates

Distal | Head width 23.5 mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80372	7	2	45 mm	23.5 mm
11-80373	7	3	53 mm	23.5 mm
11-80374	7	4	61 mm	23.5 mm
11-80375	7	5	69 mm	23.5 mm
11-80376	7	6	77 mm	23.5 mm
11-80377	7	7	85 mm	23.5 mm
11-80378	7	8	93 mm	23.5 mm
11-80379	7	9	101 mm	23.5 mm
11-80380	7	10	109 mm	23.5 mm
11-80381	7	11	117 mm	23.5 mm

Right

Radius XLOC® Poly Plates

Distal | Head width 23.5 mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80272	7	2	45 mm	23.5 mm
11-80273	7	3	53 mm	23.5 mm
11-80274	7	4	61 mm	23.5 mm
11-80275	7	5	69 mm	23.5 mm
11-80276	7	6	77 mm	23.5 mm
11-80277	7	7	85 mm	23.5 mm
11-80278	7	8	93 mm	23.5 mm
11-80279	7	9	101 mm	23.5 mm
11-80280	7	10	109 mm	23.5 mm
11-80281	7	11	117 mm	23.5 mm

Left

Radius XLOC® Poly Plates

Distal | Head width 26.8mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80572	7	2	45mm	26.8mm
11-80573	7	3	53mm	26.8mm
11-80574	7	4	61mm	26.8mm
11-80575	7	5	69mm	26.8mm
11-80576	7	6	77mm	26.8mm
11-80577	7	7	85mm	26.8mm
11-80578	7	8	93mm	26.8mm
11-80579	7	9	101mm	26.8mm
11-80580	7	10	109mm	26.8mm
11-80581	7	11	117mm	26.8mm

Right

Radius XLOC® Poly Plates

Distal | Head width 26.8mm | Titanium



ITEM NO.	HOLES HEAD	HOLES SHAFT	LENGTH	HEAD WIDTH
11-80472	7	2	45mm	26.8mm
11-80473	7	3	53mm	26.8mm
11-80474	7	4	61mm	26.8mm
11-80475	7	5	69mm	26.8mm
11-80476	7	6	77mm	26.8mm
11-80477	7	7	85mm	26.8mm
11-80478	7	8	93mm	26.8mm
11-80479	7	9	101mm	26.8mm
11-80480	7	10	109mm	26.8mm
11-80481	7	11	117mm	26.8mm

RADIUS XLOC® POLY PLATES SET


Thanks to the integrated silicone holders, the Radius XLOC® Poly Bone Plate Set can be individually configured.



Radius XLOC® Poly Plates Tray, empty

TRAY-0600-00

RADIUS XLOC® POLY INSTRUMENT SET



ITEM NO.	ITEM	DESCRIPTION	NO.
Tray-0650-00	Radius XLOC® Poly Instrument Tray		1
1-0005	Drill Bit Ø 1.9mm AO		2
1-0343	Depth gouge for screws Ø2.7 – 4.0mm		1
1-0416	Universal drill guide Ø2.5mm / Ø3.5mm		1
1-0430	Threaded drill guide Ø1.9 XLOC®		2
1-0432	Drill guide Ø1.9mm, XLOC® Poly		1
1-0610	Screw driver T10		1
1-0756	Screw driver T10 AO		2
1-0780	Screw forceps		1
1-0848	T-Handle AO		1
1-2820	Torque limiter 0.8 Nm AO	optional	1
1-2821	Torque limiter 0.8 Nm AO	optional	1
5-0120	Guide wire Ø1.6 mm L 125mm	trocar / round	6



Drill Bit Ø 1.9 mm | AO

1-0005



Depth gauge for screws Ø2.7 – 4.0 mm

1-0343



Universal drill guide Ø2.5 mm / Ø3.5 mm

1-0416



Threaded drill guide Ø1.9 XLOC®

1-0430



Drill guide Ø1.9 mm, XLOC® Poly

1-0432



Screw driver T10

1-0610



Screw driver T10 | AO

1-0756



Screw forceps

1-0780



T-Handle | AO

1-0848



Torque limiter 0.8 Nm | AO

1-2820




Torque limiter 0.8 Nm | AO

1-2821



Guide wire Ø1.6 mm | L 125 mm | trocar / round

5-0120



SiS Medical GmbH

Max-Planck-Strasse 9/6
78549 Spaichingen, Germany

Phone +49 7424 96 99 390
Fax +49 7424 96 99 399
Mail info@sismedical.de

